

CERTIFICATE OF ANALYSIS

PRODUCT	PO5090A	CHOCOLATE AGAR + VITOX
LOT NUMBER	6221774	
EXPIRY DATE	2025.03.18	
TEST DATE	2025.01.14	
REPORTING DATE	2025.01.17	

General Characteristics	Results	Specification
Colour	Conforms	Fawn brown
Appearance	Conforms	Opaque
pH	7.1	7.0 -7.4
Packaging / Presentation	Conforms	Conforms to specification
Cont. check @ 20-25 & 30-35°C for >= 72h	Conforms	Within Limits

Microbiological Performance	Control c.f.u	Test Result	Specification
Haemophilus influenzae Typ b ATCC®10211	84	84	3-5mm, cream shiny colonies
Neisseria gonorrhoeae ATCC®49226	<10000	Conforms	Good growth, cream shiny colonies
Staphylococcus epidermidis ATCC®12228	<10000	Conforms	Good growth, white colonies

For positive strains, colony count is greater than or equal to 70% of the control medium.

The information given is believed to be correct. However, both the information and the product are offered without warranty for any specific application other than that specified. The results reported were obtained at the time of release.

This certificate is produced electronically and valid without a signature

The quality control methods meet requirements of ISO 11133.