	Document Owner Department: QC	MBD-BT-SPEC-0362
		Page 1 of 2
<b>OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION</b>		
<b>X FACTOR DISCS DD0003</b>		

**X FACTOR DISCS**

**DD0003**

**Description**


Discs impregnated with Haemin for the differentiation of the Haemophilus group of bacteria.

**Microbiological Characteristics**

Place discs on to Blood Agar Base (CM0055) plates inoculated with the test organisms.

**Reactions after incubation at 37°C for 18 hours**

<i>Haemophilus influenzae</i>	ATCC® 19418	No growth
<i>Haemophilus influenzae</i>	ATCC® 49247	No growth
<i>Haemophilus parainfluenzae</i>	ATCC® 33392	No growth
<i>Haemophilus haemoglobinophilus</i>	NCTC8540	Growth

	Document Owner Department: QC	MBD-BT-SPEC-0362
		Page 2 of 2
<b>OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION</b>		
<b>X FACTOR DISCS DD0003</b>		

### Revision History

Section / Step	Description of Change	Reason for Change	Reference
Entire document	Update to new document format and correction of minor/typographical errors	N/A	N/A
Microbiological characteristics	Removal of <i>H. influenzae</i> ATCC®9334. Addition of <i>H. haemoglobinophilus</i> NCTC8540.	Change control	MOC-2023-0085