

CERTIFICATE OF ANALYSIS

PRODUCT	PO5072A	LEGIONELLA BCYEa WITH L-CYST.
LOT NUMBER	6177884	
EXPIRY DATE	2025.03.25	
TEST DATE	2024.09.24	
REPORTING DATE	2024.09.27	

General Characteristics	Results	Specification
Colour	Conforms	Jet black - Traffic black
Appearance	Conforms	Opaque
pH	6.9	6.6 -7.0
Packaging / Presentation	Conforms	Label & Print check
Cont. check @ 20-25 & 30-35°C for >= 72h	Conforms	Within Limits
Membrane filter 3	Lot 1389819	Product Code NG147-0045

Microbiological Performance	Control c.f.u	Test Result	Specification
Strain tested by membrane filtration method with above listed lot(s)			
Legionella pneumophila ATCC®33152 WDCM 00107	56	65	2- 6mm, grey-blue colonies
Strain tested by spread plate method			
Leg.pneumophila subsp.fraseri ATCC®33156 WDCM 00180	89	86	2- 4mm, grey-blue colonies

Biological performance tested in accordance with current version of ISO 11133.

The formulation of the medium conforms to ISO11731-1/2.

For positive strains, colony count is greater than or equal to 70% of the control medium.

The information given is believed to be correct. However, both the information and the product are offered without warranty for any specific application other than that specified. The results reported were obtained at the time of release.

This certificate is produced electronically and valid without a signature

The quality control methods meet requirements of ISO 11133.



The testing laboratory of Oxoid Deutschland GmbH is accredited by the German accreditation body DAkkS according to DIN EN ISO/IEC 17025 for the performance testing of media for microbiology to DIN EN ISO11133 and registered under D-PL-20190-01-00.