

CERTIFICATE OF ANALYSIS

PRODUCT	PO5323C	VIOLET RED BILE GLUCOSE AGAR (CONTACT)
LOT NUMBER	6186544	
EXPIRY DATE	2025.01.06	
TEST DATE	2024.10.14	
REPORTING DATE	2024.10.17	

General Characteristics	Results	Specification
Colour	Conforms	Signal violet
Appearance	Conforms	Transparent
pH	7.6	7.2 -7.6
Packaging / Presentation	Conforms	Label & Print check
Cont. check @ 20-25 & 30-35°C for >= 72h	Conforms	Within Limits

Microbiological Performance	Control c.f.u	Test Result	Specification
Salmonella Typhimurium ATCC®14028 WDCM 00031	83	75	2-3mm, rose cols with bile precipitate
Escherichia coli ATCC®8739 WDCM 00012	82	72	3-7mm, violet cols with bile precipitate
Enterococcus faecalis ATCC®29212 WDCM 00087	>100,000	Conforms	Total inhibition

Biological performance tested in accordance with current version of ISO 11133.

Tested in accordance with BP/EP/JP/USP. Clearly visible growth within 3 days for bacteria.

Formulation conforms to DIN 10164 and ISO 21528-1 and -2.

For positive strains, colony count is greater than or equal to 50% of the control medium.

The information given is believed to be correct. However, both the information and the product are offered without warranty for any specific application other than that specified. The results reported were obtained at the time of release.

This certificate is produced electronically and valid without a signature

The quality control methods meet requirements of ISO 11133.



The testing laboratory of Oxoid Deutschland GmbH is accredited by the German accreditation body DAkkS according to DIN EN ISO/IEC 17025 for the performance testing of media for microbiology to DIN EN ISO11133 and registered under D-PL-20190-01-00.