

	Document Owner Department: QC	BT-SPEC-0107
		Page 1 of 3
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
LYSINE MEDIUM CM0191		

LYSINE MEDIUM

CM0191

Typical Formula*

	grams per litre
Glucose	44.5
Potassium dihydrogen phosphate	1.78
Magnesium sulphate	0.89
Calcium chloride fused	0.178
Sodium chloride	0.089
Adenine	0.00178
DL-Methionine	0.000891
L-Histidine	0.000891
DL-Tryptophan	0.000891
Boric acid	0.0000089
Zinc sulphate	0.0000356
Ammonium molybdate	0.0000178
Manganese (II) sulphate	0.0000356
Iron (II) sulphate	0.0002225
L-Lysine	1.0
Inositol	0.02
Calcium pantothenate	0.002
Thiamine hydrochloride	0.0004
Pyridoxine hydrochloride	0.0004
p-Aminobenzoic acid	0.0002
Nicotinic acid	0.0004
Riboflavin	0.0002
Biotin	0.000002
Folic acid	0.000001
Agar	17.8

* adjusted as required to meet performance standards

Directions

Suspend 6.6g in 100ml distilled water. Add 1ml Potassium Lactate 50% (SR0037). With frequent agitation, bring to the boil to dissolve completely. Cool to 50°C and aseptically add the appropriate amount** of Lactic Acid 10% (SR0021) to adjust to pH 4.8 ± 0.2. Mix well and pour into sterile Petri dishes. DO NOT AUTOCLAVE. DO NOT OVERHEAT.

** This figure varies - see CM0191 label

ThermoFisher S C I E N T I F I C	Document Owner Department: QC	BT-SPEC-0107
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
LYSINE MEDIUM CM0191		

Physical Characteristics

White, free-flowing powder
 Colour on reconstitution - off white
 Moisture level - less than or equal to 7%
 pH - 5.0 ± 0.2 at 25°C
 pH - 4.8 ± 0.2 at 25°C (after adjustment)
 Clarity - clear
 Gel strength - firm, comparable to 17.8g/litre of agar

Microbiological Tests Using Optimum Inoculum Dilution

Control Medium: Sabouraud Dextrose Agar

Tested using acidified medium

Reactions after aerobic incubation at 25°C for 4 days

Medium is challenged with 1E+04 to 1E+06 colony-forming units

<i>Saccharomyces (carlsbergensis) uvarum</i>	ATCC®2700	Inhibited as slight, background film
<i>Saccharomyces cerevisiae</i>	ATCC®9763	Inhibited as slight, background film
<i>Pichia fermentans</i>	ATCC®10651	White colonies
<i>Candida vini</i>	NCYC331	White colonies

Washed, pure and mixed cultures of the listed organisms are inoculated on to the surface of the medium.

Growth of *Pichia fermentans* and *Candida vini* shall be comparable to the standard medium after incubation at 25°C for 4 days.

	Document Owner Department: QC	BT-SPEC-0107
		Page 3 of 3
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
LYSINE MEDIUM CM0191		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
Entire Document	Update to new document format and correction of typographical/minor errors. Addition of Control Medium and Result Criteria.	Change control	BT-CC-1931
Physical Characteristics	Addition of pH - 5.0 ± 0.2.	Change control	BT-CC-1931