

**CERTIFICATE OF ANALYSIS**

<b>PRODUCT</b>	<b>PO5094C</b>	<b>SAB GLUC CHLORAMPHENICOL AGAR (CONTACT)</b>
<b>LOT NUMBER</b>	6186504	
<b>EXPIRY DATE</b>	2025.01.07	
<b>TEST DATE</b>	2024.10.15	
<b>REPORTING DATE</b>	2024.10.21	

General Characteristics	Results	Specification
Colour	Conforms	Honey yellow
Appearance	Conforms	Transparent
pH	5.7	5.4 -5.8
Packaging / Presentation	Conforms	Label & Print check
Cont. check @ 20-25 & 30-35°C for >=120h	Conforms	Within Limits

Microbiological Performance	Control c.f.u	Test Result	Specification
Candida albicans ATCC®10231	80	86	2-3mm, cream colonies
Aspergillus brasiliensis ATCC®16404	23	19	10-30mm, white mycelium, black spores
Escherichia coli ATCC®8739	>100,000	Conforms	No growth
Staphylococcus aureus ATCC®6538	>100,000	Conforms	No growth

For positive strains, colony count is greater than or equal to 50% of the control medium.

The information given is believed to be correct. However, both the information and the product are offered without warranty for any specific application other than that specified. The results reported were obtained at the time of release.

This certificate is produced electronically and valid without a signature

The quality control methods meet requirements of ISO 11133.